

## What's Hot and What's Getting Hotter

BY HOWARD S. GLAZER, DDS, FAGD

### MATERIAL

#### **iBOND® Total Etch**

Heraeus Kulzer, LLC  
300 Heraeus Way  
South Bend, IN 46614-2517  
800.431.1785  
[www.heraeusdentalusa.com](http://www.heraeusdentalusa.com)

From the company that first introduced and led the way to seventh generation adhesives with iBond® comes a great fifth generation adhesive. Heraeus Kulzer's iBOND Total Etch offers great strength, provides for desensitization, and is easy to handle. This is not a step backward but rather a giant step forward in advancing adhesive dentistry. What makes this fifth generation product unique is the nano-filler composition, proprietary to Heraeus, that provides for

greater bond strength to the enamel and dentin. iBOND Total Etch is a two-step bonding system, and application is simple: First, etch and rinse; then, apply one single layer of iBOND Total Etch. No agitation is required. The material air-dries quickly and will not be blown away by an air stream because it has a thicker consistency than most other similar products. No refrigeration is required, and it comes in unit-dose packaging as well as in a multi-dose 4 ml bottle. As with most adhesive agents, iBOND Total Etch is indicated for bonding all indirect and direct composite restorations and can treat hypersensitive areas of teeth such as cervical erosion. This product also can be used for bonding porcelain laminate veneers. iBOND products continue to be innovative in providing great bond strengths and easier application, which benefits both the doctor and the patient.




## iBOND®

Proven to outperform the market leader. Every time.

- High bond strength to both dentin and enamel
- Clinically proven to provide reliable and consistent results
- Convenient handling with easy application
- Demonstrates reduction in postoperative sensitivity

**iBOND® Total Etch is clinically proven to outperform OptiBond® Solo Plus in bond strength.**

In a study conducted by REALITY, iBOND Total Etch was shown to have higher bond strengths to dentin and enamel than OptiBond Solo Plus.



Source: REALITY, October 2009. Sponsored by Heraeus. Data on file. OptiBond Solo Plus is not a registered trademark of Heraeus Kulzer GmbH.

Average standard of deviation reveals iBOND Total Etch to have a higher, more consistent bond strength to both enamel and dentin.

**Purchase 2  
iBOND® Self Etch  
or iBOND® Total  
Etch, receive 1 of  
the same FREE!**  
Call 1-877-258-5937  
for details.



iBOND®  
Bonding with Confidence